

EXPRE^S²ION

BIOTECHNOLOGIES

EXPRES²ION BIOTECHNOLOGIES

- Founded in 2010 in Copenhagen, Denmark
- A CRO spin-off from Affitech A/S
- Proprietary technology platform for production of proteins based on *Drosophila* S2 insect cells
- > 10 years experience with S2 cells from research to GMP production for human clinical trials

OUR FOCUS


- 100% focus on Drosophila S2 cells
- Vector and cell line development
- Upstream process development
- Small-scale Protein production
- Process Transfer to CMOs for GMP production

AIM

- Introduce the *Drosophila* S2 system as a complement to the Baculovirus system
- Make the S2 system available to academic centers for free



A COMPLETE PROTEIN EXPRESSION SYSTEM

- >10 years experience with *Drosophila* S2
- Vaccine clinical trials up to Phase II



DROSOPHILA S2 INSECT CELLS

- Derived from the late stages of *Drosophila melanogaster* (fruit fly) embryos (Schneider, 1972)
- Doubling times of 18-24 hrs (both serum-containing and serum-free medium)
- Cell size: 9-10 μm in diameter
- Grow well between 5 and 140E6 cells/ml in suspension
- Grow at 23°C to 27°C
- Very robust



LITERATURE EXAMPLES OF S2 EXPRESSED PROTEINS

	Family	S2 Produced
Secreted	Anti-angiogenic	Endostatin, Canstatin
	Cytokines & growth factors	VEGF, IL5, IL12, EPO, RANKL, TNF, lymphotactin
	Enzymes	Sea raven type II antifreeze protein
	Plasma proteins	Human transferrin
	Monoclonal antibodies	IgG, IgM, scFv, Fc, chimeras & Fab fragments
Viral	Viral Proteins	JEV, rabies virus, West Nile virus, Influenza, Dengue fever
	Virus like particles	JEV glycoprot (E), HIV, Rotavirus
Complex surface	Receptors	CD23*, EGFR*, Her2*, ADAM33*, glucagon rec., MHC class II I-E _d
	Integrins & related receptors	α IIb β 3, suPAR
	GPCRs	D1 and D2 dopamine receptors, human mu opioid receptor, IL-5 receptor α , chemokine CxCR3, CXCL11 receptors

*soluble versions

EXPRES²ION

FURTHER EXAMPLES OF S2 EXPRESSED PROTEINS

Protein class	Protein name	Source	Reference	Comments
GPCR ^a	M ₁ muscarinic acetylcholine receptor	<i>Drosophila</i>	Millar et al. (1995) Yagodin et al. (1999)	Robust receptor expression; ~2.5 pmol/mg Imaging of intracellular Ca ²⁺
GPCR	M ₂ muscarinic receptor	Pig	Graziano et al. (1998)	Robust receptor expression; ~5 pmol/mg
GPCR	D _{2L} dopamine receptor	Rat	Schetz et al. (2003)	Robust receptor expression; ~10 pmol/mg. Stereoselective drug binding is similar to that in mammalian cells.
GPCR	D ₁ dopamine receptor	Rat	Schetz et al. (2003)	Robust receptor expression; ~8.5 pmol/mg. Stereoselective drug binding is similar to that in mammalian cells.
GPCR	Glucagon receptor	Human	Graziano et al. (1998)	Ultra-high expression; ~100 pmol/mg. Anisotropy of fluorescent ligands.
GPCR	Tachykinin receptor (STKR)	<i>Stomoxys calcitrans</i>	Torfs et al. (2000) Torfs et al. (2001)	Positively coupled to phospholipase C and adenylyl cyclase Receptor expression; ~0.15 pmol/million cells. Imaging of intracellular Ca ²⁺
GPCR	Octopamine receptor (OMAB)	<i>Drosophila</i>	Van Poyer et al. (2001)	Positively coupled to adenylyl cyclase and activated by the phenolamines: tyramine and octopamine
GPCR	Leucokinin receptor	<i>Drosophila</i>	Radford et al. (2002)	Immunocytochemistry and western blotting
GPCR	Inositol (1,4,5)-triphosphate receptor (IP ₃ R)	<i>Drosophila</i>	Pollock et al. (2003)	Ultra-high protein expression; ~400 mg/liter
GPCR	Interleukin 5 (IL-5) receptor	Human	Johanson et al. (1995)	Robust receptor expression; 17 mg/liter. The peptide ligand, IL-5, is also expressed at ~22 mg/liter.
GPCR	Chemokine(CXCR3) ^b	Human	Hensbergen et al. (2001)	Expression of CXCR3 receptors is similar to that in mammalian cells
GPCR	Chemokine(CXCL11) ^f	Human/ murine	Hensbergen et al. (2001)	Expression of CXCL11 receptors is similar to that in mammalian cells
Ion channel	γ-Aminobutyric acid receptor (GABA _A)	<i>Drosophila</i>	Millar et al. (1994)	Robust expression; ~3 pmol/mg Whole cell patch clamp Single-channel recordings.
Ion channel	γ-Aminobutyric acid receptor (GABA)	<i>Drosophila</i>	Grolleau et al. (2000)	
Ion channel	Nicotinic acetylcholine receptor	<i>Drosophila</i>	Lansdell et al. (1997)	Robust receptor expression; ~6.5 pmol/mg
Ion channel	Chloride channel	<i>Drosophila</i>	Asmild et al. (2000)	Endogenous channels. Whole cell patch clamp recording gives an estimate of 0.02 channels/μm ² . Patch clamp recordings
Ion channel	Voltage-dependent sodium channel	<i>Drosophila</i>	Meadows et al. (2002)	
Ion channel	Transient receptor potential-like (TRPL)	<i>Drosophila</i>	Goel et al. (2001)	Immunoprecipitation and immunoblotting. Imaging of intracellular Ca ²⁺ . Comparison of expression with <i>S9</i> cells.
Enzyme	Chloramphenicol 3-O-acetyltransferase	Bacterial	Di Nocera and Dawid (1983)	Effective reporter protein
Enzyme	Galactokinase	Bacterial	Johansen et al. (1989)	Robust expression; ~3 mg/liter. Effective reporter protein.
Enzyme	Dihydrofolate reductase	Bacterial	Van der Straten et al. (1987)	Confers resistance to methotrexate
Enzyme	Hygromycin B phosphotransferase	Bacterial	Van der Straten et al. (1987)	Confers resistance to hygromycin B
Enzyme	Neomycin	Bacterial	Van der Straten et al. (1987)	Confers resistance to G418
Enzyme	Glucosidase II	<i>Drosophila</i>	Ziak et al. (2001)	Endogenous enzyme and only one isoform expressed

Protein class	Protein name	Source	Reference	Comments
Peptide	Neurotensin	Human	Friry et al. (2002)	Levels of mRNA transcribed similar to that in mammalian cells
Antibody	IgG ₁ monoclonal antibody	Human	Kirkpatrick et al. (1995a)	BiP-mediated assembly of both heavy and light chains
Oncogene protein	H-ras	Human	Johansen et al. (1989)	H-ras is growth-inhibitory, thus an inducible promoter is required
Cytoskeletal protein	Cadherin	<i>Drosophila</i>	Dubreuil et al. (1999)	Western blotting
Cytoskeletal protein	Neuroglian	<i>Drosophila</i>	Dubreuil et al. (1999)	Western blotting
Glycoprotein	HIV-gp120 ^d	AIDS virus	Culp et al. (1991)	Robust protein expression; ~35 mg/liter. Glycosylation is limited to simple core sugars.
Glycoprotein	Erythropoietin (EPO)	<i>Drosophila</i>	Lee et al. (2000)	Robust protein expression is; ~12mg/liter
Apoptotic protein	Presenilin	<i>Drosophila</i>	Ye and Fortini (1999)	Western blotting
Apoptotic protein	Inhibitor of apoptosis (IAP1)	<i>Drosophila</i>	Muro et al. (2002)	Comparison of expression with that in <i>S9</i> cells

^aGPCR: G protein-coupled receptor.
^bThe alternative nomenclature for the chemokine (CXCR3) receptor is interferon inducible protein-10 (IP-10).
^cThe alternative nomenclature for the chemokine (CXCL11) receptor is interferon inducible protein-9 (IP-9).
^dHIV: Human immunodeficiency virus.

*“Following drug selection, a polyclonal S2 cell line can be induced to express on the order of 2 to 100 pmol/mg membrane protein for G-coupled protein receptors, 4000 to 100,000 sites/cell for other membrane receptors and 3 to 35 mg/liter for soluble and secreted proteins.” **

*** Protein Expression in the *Drosophila* Schneider 2 Cell System, Current Protocols in Neurosciences, John A. Schetz, Eswar P.N. Shankar, 2004**


Proteins expressed by ExpreS²ion

Protein Class

Secreted Binding Proteins
Biomarkers
Cytokines
Differentiation Factors
Growth factors
HDACs
IgGs
Kinases
Malaria antigens
Receptors *
Membrane proteins *
Metal binding proteins
Methyltransferases
Plasma Proteins
Proteases
Secreted Signaling Molecules
Viral antigens


* Soluble version

Proteins expressed using ExpreS²




IMPROVED YIELDS

Vector development




Medium development



*Adaptation to individual media were performed

QUICK TO PROTEIN AND STABLE PROCESS

Transfection



EXPRES²ION

BENEFITS of S2 SYSTEM

Compared to Baculovirus insect cell systems


- Non-viral, Non-lytic
- More homogeneous product
- High reproducibility and consistency between batches
- Potential for Higher yield
- Long-term stability of frozen cell stocks

Compared to other S2 systems

- Increased yield
- Protocols for large-scale production
- Including Perfusion processes

EXPRES²ION

More homogeneous product




BEVS: Baculovirus Expression Vector System

EXPRES²ION

THE EASY WAY TO PROTEINS

High-throughput protein expression in *Drosophila* S2 cells

- Robust cells and increased expression allows growth, transfection and expression testing in 1ml, 96 well format
- Fast scale-up for selected proteins for structure and functional studies



Take-home messages

- Complementary system to Baculovirus for structure and functional studies
- Opportunity to further improve hit-rate for complex proteins
- Potential yield and consistency advantages for secreted proteins


Questions?

GLYCOSYLATION

Insect (S2)

CHO

Human



EXPRES²ION

N-glycosylations found in RANKL based on MALDI-TOF analysis

0 N-glycosylation

+1 N-glycosylation

+2 N-Glycosylation

